









STEPS
Sunrise Training & Education Programs

Power Technical School

Presented by:
Ken Kalinowski, Senior Service Technician
Sunrise Medical LLC

MANUALPOWERADULTPEDIATRICSEATINGGERIATRICSCONTROLSFUNDING

Agenda and Course Objectives

8:00 Power Wheelchair Details and Overview
Examples: Pulse, S636, QM-710
Batteries & Access
Suspension – Electronics
ASAP Seating Adjustments


9:30 Electronics
Battery Chargers
VR2 – R-Net 90

10:00 Break Time

10:15 Programming
VR2 with Qtronix Pad; R-Net
Programming – DTT & On Board Programming

11:45 Discussion and Overview

12:15 End of Half Day




CONFIDENTIAL AND PROPRIETARY
Duplication or Distribution Prohibited
6/20/12


2

Rehab Power Ladder Example


- For Mid Wheel Drive (MWD)




Pulse 6 –
BC – K0848




QM-710 w/
Rnet BC – K0848




Pulse 6 -
SC – K0856




QM-7 Series Rnet
SC – K0856



QM-710 w/ Rnet
MPC – K0861





CONFIDENTIAL AND PROPRIETARY
Duplication or Distribution Prohibited
6/20/12

3

Mid & Rear Wheel Drive Power Family



Pulse 6
MWD



QM-710
MWD



S-636
RWD

All Group Three Power Chairs




CONFIDENTIAL AND PROPRIETARY
Duplication or Distribution Prohibited
6/20/12

4

Example: Power Portfolio Strategy


- Simple and reliable
 - Focused
 - Example: Specializing in Group 3 Power with three models
 - Proven Reliability in Components
 - Motors
 - Electronics
 - Tilt Module
 - Made in the USA




CONFIDENTIAL AND PROPRIETARY
Duplication or Distribution Prohibited
6/20/12

5


Example: Models



Pulse 6 BC
K0848



Pulse 6 SC
K0856



CONFIDENTIAL AND PROPRIETARY
Duplication or Distribution Prohibited
6/20/12

6

MWD - Pulse 6 Specifications

- HCPCS Code
 - Pulse 6 BC: K0848
 - Pulse 6 CC: K0849
 - Pulse 6 SC: K0856
- Drive Wheel Position
 - MWD
- Speed
 - 6.0 mph
- Drive Wheel Tires
 - 13"
- Battery Type
 - 22 NF
- Turning Radius
 - 22"
- Product Weight
 - 235 - 260 lbs.
- User Weight Capacity
 - 300 lbs.
- Overall Base Width
 - 24.5"
- Overall Length
 - 33.9"
- Seat Width
 - 12" – 22"
- Seat Depth
 - 12" – 22"
- Seat to Floor Height
 - BC: 16.5", 17.5", 18.5"
 - CC: 17", 18", 19"
 - SC: 17.5", 18", 19.5"



CONFIDENTIAL AND PROPRIETARY
Duplication or Distribution Prohibited
6/20/12
7

VR-2 Electronics

- Highly reliable and fully functioning non-expandable control platform
- Economical thru-drive capability for tilt
- Attendant control option



CONFIDENTIAL AND PROPRIETARY
Duplication or Distribution Prohibited
6/20/12
8

Two Point Reclining Armrest

- Provides additional support during transfers
- Ability to unlatch to flip back
- Can be adjusted to orient the arm parallel to the ground
- Remains solid during tilt
- Acts as an additional thigh support




CONFIDENTIAL AND PROPRIETARY
Duplication or Distribution Prohibited
6/20/12
9

Pulse 6 : Footprint


Performance

- Compact Footprint



24"

33.9"

STEPS 
Senior Training & Education Program


CONFIDENTIAL AND PROPRIETARY
Duplication or Distribution Prohibited
6/20/12

10


Pulse 6 : Performance

Performance

- 6 mph standard
- Dynamic Suspension



18°


STEPS 
Senior Training & Education Program

CONFIDENTIAL AND PROPRIETARY
Duplication or Distribution Prohibited
6/20/12

11


Pulse 6: Backrest Options

- **JAY SPO & ART SPO Backrests**
 - Great complement to ASAP and SPOT adjustable seat frames
 - Modular foam insert construction
 - Countless in-field adjustment opportunities



JAY SPO Back

A.R.T. SPO Back

STEPS 
Senior Training & Education Program

CONFIDENTIAL AND PROPRIETARY
Duplication or Distribution Prohibited
6/20/12

12

Pulse 6 : Tilt

Single Power Option Tilt (SPOT)

- As low as 17" STFH
- 50° range



Pulse 6: Hanger Options

- Fixed Centermount
- 90 deg
- 75 deg taper
- 70 deg Swing-away
- 65 deg Swing-away
- Power Center mount single and Ped.



ASAP II Seating Adjustments

- The ASAP II seating system is fully adjustable to allow changes from 16" to 20" or
- 18" – 22" (requires -Frame Weldment and Seat Pan Plates change) - for seat growth.
- The picture shown in figure 1 – shows a 20" wide seat width
This detail will show how easy to change back to 16" wide.



ASAP II Seating Adjustments

- **Step 1** – Confirm that you have a 3mm – 4mm – and 5mm Allen Hex keys available for making adjustments.
- **Step 2** – Loosen front M8 Cap screw bolt – with a 4mm Allen Button Head Cap screw driver or Allen key. These are seen attached to plate on Figure 1. There is no need to take both out for this plate slides with size changes of seat.
- **Step 3** – Remove 3mm Allen Hex Key Socket Head Cap Screws M8 x 1.25 size that are supporting the wings left and right to the seat frame



CONFIDENTIAL AND PROPRIETARY
Duplication or Distribution Prohibited
6/20/12
16

ASAP II Seating Adjustments

- **Step 4** – Take out the 5mm head bolts that are used to hold bottom seat in place – totally remove out of its position. Shown on figure below.
- **Step 5** – Take out bolts that are in Rear seat weldment – brace using a 4mm Allen Hex Key and a 10 mm open end wrench or ratchet shown below in pictures.



CONFIDENTIAL AND PROPRIETARY
Duplication or Distribution Prohibited
6/20/12
17

ASAP II Seating Adjustments

- **Step 6** – Loosen Back bracket off of chair to allow for changes to seat – note it is recommended to use the back width to accommodate the user – so 16 with 16 wide – 18 with 18" and so fourth.
Note: However if the need occurs that a 1/2" change is needed the brackets shown in picture in Figure - show you have the 1/2" change needed for a 1" change .
- **Step 7** – Remove 2 bolts from Seat Bracket attached to the back cane assembly – with a 5mm Allen Hex Key and 1/2 inch nut using either open end or Ratchet.



CONFIDENTIAL AND PROPRIETARY
Duplication or Distribution Prohibited
6/20/12
18

ASAP II Seating Adjustments

- Push Frame pieces back into place
- Tighten 5mm bolts of Seat Base – tighten 4mm bolts of seat base only one side on front and back this plate will slide. Also tighten the bolts from each of the wings where it mounts to seat -using 3mm.
- Finally Tighten the 3mm Allen Hex Key Button Head Screws back into place



CONFIDENTIAL AND PROPRIETARY
Duplication or Distribution Prohibited
6/20/12
19

ASAP II Seating Adjustments

- Retighten Front and Back 4mm Allen Hex Key Socket Head Cap Screws in place for support piece.
- Retighten in place the bolts supporting the back.



CONFIDENTIAL AND PROPRIETARY
Duplication or Distribution Prohibited
6/20/12
20

ASAP II Seating Adjustments

- Check all bolts and if you run into any of the bolts not matching see bottom of Seat Frame depicted below 2 adjustment area of inside where frame inserts and outside where frame attaches.



CONFIDENTIAL AND PROPRIETARY
Duplication or Distribution Prohibited
6/20/12
21

ASAP II Seating Adjustments

- Check to make sure all parts are tightened and check width 16" wide complete and measures up to expectations




Pulse 6: Tilt Module

- Tilt Module mounts between the ASAP seat and power base
- The module is replaced as an entire unit
- Tilt modules can be added after initial purchase if necessary



Pulse 6: Tilt Module Retro Kit



- Note: If adding tilt to a chair w/o tilt, Retro Kit must be ordered and the following electronics must be ordered separately:
 - VR2 w/ thru drive (Joystick-105376 & Controller – 104451)
 - R-net or VR2 w/ Toggle (Toggle retro kit – 110842)
 - R-net w/ thru drive (no other items needed)

P/N	Item Number	Description	Comments
1	105376	VR2 w/ thru drive (Joystick-105376 & Controller – 104451)	
2	105376	VR2 w/ thru drive (Joystick-105376 & Controller – 104451)	
3	105376	VR2 w/ thru drive (Joystick-105376 & Controller – 104451)	
4	105376	VR2 w/ thru drive (Joystick-105376 & Controller – 104451)	
5	105376	VR2 w/ thru drive (Joystick-105376 & Controller – 104451)	

Mounting Tilt to ASAP Seat

FRONT

FRONT TILT INTERFACE PLATE

FRONT TILT INTERFACE BRACKET

REAR BRACKET MOUNTING

REAR TILT INTERFACE BRACKET

SEAT DEPTH	FRONT PLATE MOUNTING	REAR BRACKET MOUNTING
18-20"	E-31 Y-32	K-32 Y-33
18-20"	E-31 Y-32	K-33 Y-34

STEPS
Senior Training & Education Program

CONFIDENTIAL AND PROPRIETARY
Duplication or Distribution Prohibited
6/20/12

25

Harness Routing Instructions

- For chairs with tilt install the adaptor harness.
 - The 3-way connector plugs into the controller OBC port.
 - The 6-way connector accepts the actuator harness from the tilt.
 - The 4-way connector mates either to the TM40 flying lead or the drive-through harness.

VR2-90 Controller Connector View

See the photos in upcoming slides

STEPS
Senior Training & Education Program

CONFIDENTIAL AND PROPRIETARY
Duplication or Distribution Prohibited
6/20/12

26

Serial Number Location

Old Placement of Serial Number


New Placement of Serial number

STEPS
Senior Training & Education Program


CONFIDENTIAL AND PROPRIETARY
Duplication or Distribution Prohibited
6/20/12

27

22NR Batteries and Connections




Pulse 6: Shroud Connections



*New Shroud Connections,
May 2009*

Pin	Item Number	Description	Remarks
1	104400	HARDWARE SHROUD	Includes 3 Hardware
1,2	104410	SHROUD TOP KIT RED	WEDGEABLES AND FASTENERS
1,2	104413	SHROUD TOP KIT BLACK	WEDGEABLES AND FASTENERS
1,3	104409	SHROUD FRONT KIT RED	WEDGEABLES AND FASTENERS
1,3	104411	SHROUD FRONT KIT BLACK	WEDGEABLES AND FASTENERS
1,4	104408	SHROUD REAR KIT RED	WEDGEABLES AND FASTENERS
1,4	104412	SHROUD REAR KIT BLACK	WEDGEABLES AND FASTENERS
5	112400	THUMB SCREW 1/8" x 3/4"	

Wire Routing for Main Harness



Wire Routing for Main Harness

- **Inline fuse requirement**
 - In order to satisfy a strict reading of ANSI/RESNA section 14 and properly protect the power wiring against over-current faults, a protection device must be placed as close to the battery terminal as possible. The circuit breakers we use do not lend themselves to this, so an inline fuse near the battery post is the preferred solution. These fuses are required whether the breaker is retained or not.
- **Advanced electronic controls**
 - The motor controllers that have become industry standard over the last 20 years now handle the motor overload protection function previously handled by the circuit breaker. The onboard microprocessor tracks the motor currents over time and reduces them according to parameters determined by motor and chair testing. In fact, because of PWM technology efficiency, the circuit breaker no longer sees enough current in a stall condition to properly protect the motors. This must be done by monitoring not the battery current but the actual motor currents as the controller does.



CONFIDENTIAL AND PROPRIETARY
Duplication or Distribution Prohibited
6/25/12
31

Removal of Motor Assembly



Remove bolt from front arm first



Bolt with spacer & outside washer




Bolt with spacer & inside washer
Don't reverse linkage small offset



CONFIDENTIAL AND PROPRIETARY
Duplication or Distribution Prohibited
6/25/12
32

Pulse 6: Motor Assembly



Part	Item Number	Description	Remarks
1,2,3A,10,14	10422	MOTOR ASSEMB ELECTROCRUFT L/R2	LEFT - PULSE 6 W/IR2
1,2,3A,10,14	10423	MOTOR ASSEMB ELECTROCRUFT R/R2	RIGHT - PULSE 6 W/IR2
1,2,3B,10,14	10429	MOTOR ASSEMB ELECTROCRUFT L/RW	LEFT - PULSE 6 W/IRW
1,2,3B,10,14	10430	MOTOR ASSEMB ELECTROCRUFT R/RW	RIGHT - PULSE 6 W/IRW
5	10903	MOUNT MOTOR W/ELMENT L	LEFT
5	10904	MOUNT MOTOR W/ELMENT R	RIGHT
10	10906	COO RELEASE WHEEL AND PIN	RIGHT
1,2,3A,4,14	10549	MOTOR ASSEMB L/R2 W/COO	LEFT - PULSE 6 W/IR2
1,2,3A,4,14	10549	MOTOR ASSEMB L/R2 W/COO	RIGHT - PULSE 6 W/IR2
1,2,3B,4,14	10549	MOTOR ASSEMB L/R2 W/COO	LEFT - PULSE 6 W/IRW
1,2,3B,4,14	10549	MOTOR ASSEMB L/R2 W/COO	RIGHT - PULSE 6 W/IRW
10	10907	BRUSHLESS MOTOR (10.5V)	2 Motors, 1 each




CONFIDENTIAL AND PROPRIETARY
Duplication or Distribution Prohibited
6/25/12
33

Pulse 6: Motor Assembly*

*After S/N PLS-101619

P/N	Qty	Description	Remarks
1	1	WHEEL ASSEMBLY	
2	1	WHEEL MOUNTING BRACKET	
3	1	WHEEL MOUNTING BRACKET	
4	1	WHEEL MOUNTING BRACKET	
5	1	WHEEL MOUNTING BRACKET	
6	1	WHEEL MOUNTING BRACKET	
7	1	WHEEL MOUNTING BRACKET	
8	1	WHEEL MOUNTING BRACKET	
9	1	WHEEL MOUNTING BRACKET	
10	1	WHEEL MOUNTING BRACKET	
11	1	WHEEL MOUNTING BRACKET	
12	1	WHEEL MOUNTING BRACKET	
13	1	WHEEL MOUNTING BRACKET	
14	1	WHEEL MOUNTING BRACKET	
15	1	WHEEL MOUNTING BRACKET	
16	1	WHEEL MOUNTING BRACKET	
17	1	WHEEL MOUNTING BRACKET	
18	1	WHEEL MOUNTING BRACKET	
19	1	WHEEL MOUNTING BRACKET	
20	1	WHEEL MOUNTING BRACKET	
21	1	WHEEL MOUNTING BRACKET	
22	1	WHEEL MOUNTING BRACKET	
23	1	WHEEL MOUNTING BRACKET	
24	1	WHEEL MOUNTING BRACKET	
25	1	WHEEL MOUNTING BRACKET	
26	1	WHEEL MOUNTING BRACKET	
27	1	WHEEL MOUNTING BRACKET	
28	1	WHEEL MOUNTING BRACKET	
29	1	WHEEL MOUNTING BRACKET	
30	1	WHEEL MOUNTING BRACKET	
31	1	WHEEL MOUNTING BRACKET	
32	1	WHEEL MOUNTING BRACKET	
33	1	WHEEL MOUNTING BRACKET	
34	1	WHEEL MOUNTING BRACKET	
35	1	WHEEL MOUNTING BRACKET	
36	1	WHEEL MOUNTING BRACKET	
37	1	WHEEL MOUNTING BRACKET	
38	1	WHEEL MOUNTING BRACKET	
39	1	WHEEL MOUNTING BRACKET	
40	1	WHEEL MOUNTING BRACKET	
41	1	WHEEL MOUNTING BRACKET	
42	1	WHEEL MOUNTING BRACKET	
43	1	WHEEL MOUNTING BRACKET	
44	1	WHEEL MOUNTING BRACKET	
45	1	WHEEL MOUNTING BRACKET	
46	1	WHEEL MOUNTING BRACKET	
47	1	WHEEL MOUNTING BRACKET	
48	1	WHEEL MOUNTING BRACKET	
49	1	WHEEL MOUNTING BRACKET	
50	1	WHEEL MOUNTING BRACKET	
51	1	WHEEL MOUNTING BRACKET	
52	1	WHEEL MOUNTING BRACKET	
53	1	WHEEL MOUNTING BRACKET	
54	1	WHEEL MOUNTING BRACKET	
55	1	WHEEL MOUNTING BRACKET	
56	1	WHEEL MOUNTING BRACKET	
57	1	WHEEL MOUNTING BRACKET	
58	1	WHEEL MOUNTING BRACKET	
59	1	WHEEL MOUNTING BRACKET	
60	1	WHEEL MOUNTING BRACKET	
61	1	WHEEL MOUNTING BRACKET	
62	1	WHEEL MOUNTING BRACKET	
63	1	WHEEL MOUNTING BRACKET	
64	1	WHEEL MOUNTING BRACKET	
65	1	WHEEL MOUNTING BRACKET	
66	1	WHEEL MOUNTING BRACKET	
67	1	WHEEL MOUNTING BRACKET	
68	1	WHEEL MOUNTING BRACKET	
69	1	WHEEL MOUNTING BRACKET	
70	1	WHEEL MOUNTING BRACKET	
71	1	WHEEL MOUNTING BRACKET	
72	1	WHEEL MOUNTING BRACKET	
73	1	WHEEL MOUNTING BRACKET	
74	1	WHEEL MOUNTING BRACKET	
75	1	WHEEL MOUNTING BRACKET	
76	1	WHEEL MOUNTING BRACKET	
77	1	WHEEL MOUNTING BRACKET	
78	1	WHEEL MOUNTING BRACKET	
79	1	WHEEL MOUNTING BRACKET	
80	1	WHEEL MOUNTING BRACKET	
81	1	WHEEL MOUNTING BRACKET	
82	1	WHEEL MOUNTING BRACKET	
83	1	WHEEL MOUNTING BRACKET	
84	1	WHEEL MOUNTING BRACKET	
85	1	WHEEL MOUNTING BRACKET	
86	1	WHEEL MOUNTING BRACKET	
87	1	WHEEL MOUNTING BRACKET	
88	1	WHEEL MOUNTING BRACKET	
89	1	WHEEL MOUNTING BRACKET	
90	1	WHEEL MOUNTING BRACKET	
91	1	WHEEL MOUNTING BRACKET	
92	1	WHEEL MOUNTING BRACKET	
93	1	WHEEL MOUNTING BRACKET	
94	1	WHEEL MOUNTING BRACKET	
95	1	WHEEL MOUNTING BRACKET	
96	1	WHEEL MOUNTING BRACKET	
97	1	WHEEL MOUNTING BRACKET	
98	1	WHEEL MOUNTING BRACKET	
99	1	WHEEL MOUNTING BRACKET	
100	1	WHEEL MOUNTING BRACKET	



CONFIDENTIAL AND PROPRIETARY
Duplication or Distribution Prohibited
6/20/12

34

Pulse 6: Suspension



Adjustable Shock
Front Arm Linkage
Damper



CONFIDENTIAL AND PROPRIETARY
Duplication or Distribution Prohibited
6/20/12

35

Pulse 6: Disassembly of Leaf Spring*

*Discontinued 5/09



Remove 2 outside bolts



Access to Leaf Spring




Remove 2 inside bolts & fixed bracket



CONFIDENTIAL AND PROPRIETARY
Duplication or Distribution Prohibited
6/20/12

36

Pulse 6: Suspension



Part	Part Number	Description	Remarks
1	100001	LEFT DRIVE MOTOR (VARIABLE)	LEFT
2	100002	RIGHT DRIVE MOTOR (VARIABLE)	RIGHT
3	100003	CASTER AND REAR WHEEL ASSEMBLY	REAR
4	100004	CASTER AND REAR WHEEL ASSEMBLY	LEFT
5	100005	WHEEL AND TIRE ASSEMBLY	FRONT
6	100006	WHEEL AND TIRE ASSEMBLY (12")	LEFT
7	100007	WHEEL AND TIRE ASSEMBLY (12")	FRONT
8	100008	FRONT TO REAR LINK	FRONT
9	100009	CASTER AND REAR WHEEL ASSEMBLY	REAR
10	100010	CASTER AND REAR WHEEL ASSEMBLY	LEFT
11	100011	FRONT DRIVE MOTOR	FRONT
12	100012	FRONT DRIVE MOTOR	FRONT

STEPS Senior Training & Education Programs

CONFIDENTIAL AND PROPRIETARY
Duplication or Distribution Prohibited
6/20/12

37

Pulse 6: Controller Mount



- Mounts with positive lock
- Easy removal without tools
- Locks tight into frame

STEPS Senior Training & Education Programs

CONFIDENTIAL AND PROPRIETARY
Duplication or Distribution Prohibited
6/20/12


38

Group 3 Products

K0848

Frame seating, No Tilt

Pulse 6 BC



Quickie Rhythm w/ ASAP seating


Inmacare TDx-SP w/ ASBA Seating

Pride Quantum 600 w/ Synergy Seating

K0856

Single power option

Pulse 6 SC



Quickie Rhythm w/ SPOT tilt

Inmacare TDx-SP w/ tilt

Pride Quantum 6000 w/ tilt


STEPS Senior Training & Education Programs

CONFIDENTIAL AND PROPRIETARY
Duplication or Distribution Prohibited
6/20/12


39

Rehab Power Ladder Example


Rear Wheel Drive




Quickie P222 SE
K0877




Quickie S636
K0856 & K0861



Quickie 646
K0877 & K0884





CONFIDENTIAL AND PROPRIETARY
Duplication or Distribution Prohibited
6/20/12

40

Quickie S-646 SE & P-222 SE

Updated - May 21, 2012







CONFIDENTIAL AND PROPRIETARY
Duplication or Distribution Prohibited
6/20/12

41

Changes to S646 and P222SE

- Standardize Motors and Electronics across the entire power portfolio.
- Simplify overall offering.
- Increase overall reliability.
- Update electronics to a more robust solution.
- Add an LED package so Europe can access these same products. The LED package is also becoming a popular group 4 option in the US.





CONFIDENTIAL AND PROPRIETARY
Duplication or Distribution Prohibited
6/20/12

42

S-636: Short Wheelbase







CONFIDENTIAL AND PROPRIETARY
Duplication or Distribution Prohibited
6/20/17

46


S-636: Short Wheelbase




- During tilt and recline, the dynamic rear anti-tip suspension is locked down (release levers open). Because the anti-tip is locked, the S-636 is able to support a shorter wheel base. Bottom line: shorter wheel base = greater maneuverability



Seat Down – Dynamic Anti-Tip Release Levers are depressed and anti-tip suspension is enabled




Seat Up – Dynamic Anti-Tip Release Levers are open and anti-tip suspension is disabled



CONFIDENTIAL AND PROPRIETARY
Duplication or Distribution Prohibited
6/20/17

47


S-636: RWD/MWD Hybrid



- Mid-wheel drive chairs have optimal maneuverability in tight spaces

The S-636 is a rear biased mid-wheel hybrid that offers the best of both worlds

- Rear wheel drive chairs track straight at high speeds and have good curb climbing ability



CONFIDENTIAL AND PROPRIETARY
Duplication or Distribution Prohibited
6/20/17

48

S-636: Rear Tail Gate

- The S-636 has one of the most easily accessible battery compartments in the power chair market. It begins with:
 - A toggle clamp lock – easy to open, positive pressure lock
 - Fold down rear gate – supports battery tray



CONFIDENTIAL AND PROPRIETARY
Duplication or Distribution Prohibited
6/25/12
49

S-636: Servicing Batteries

- Batteries can be accessed without the user having to leave the chair.



CONFIDENTIAL AND PROPRIETARY
Duplication or Distribution Prohibited
6/25/12
50

S-636: What's New or Changing

- S-636 changed to R-net Electronics in May 2010
- Discontinuation of Smart Seat tilt and recline system – will keep tilt and move to current recline platform used on the Quickie Rhythm.
- Moved to Linux motors





CONFIDENTIAL AND PROPRIETARY
Duplication or Distribution Prohibited
6/25/12
51

S-636: Old Motors

Lever Arm Release

Old Motors

Gear Release Freewheel



Changing: Difficult to engage, often requires rocking to release

Changing: Complexity of gear box makes service more difficult. Also, occasionally leaks grease.

STEPS

CONFIDENTIAL AND PROPRIETARY
Duplication or Distribution Prohibited
6/25/12
52

S-636: New Motors


- Communication
 - No change in performance or price.
 - Only noticeable change is in the freewheel method
- New motor continues our commitment to make the product simple and reliable

STEPS

CONFIDENTIAL AND PROPRIETARY
Duplication or Distribution Prohibited
6/25/12
53

S-636: New Motors

Cog Release – Old Style	Brake Release – New Style
Push Force easier	Push Force greater
Difficult to engage	Easy to engage
Hand Operated only	Hand or Foot operated
Prone to leaking	Gearbox completely sealed
More moving parts	Less parts.
More complex	Simple




STEPS


CONFIDENTIAL AND PROPRIETARY
Duplication or Distribution Prohibited
6/25/12
54

S-636: New Motors


- Linux 4 pole motors

Part Number	Description
107247	Motor 10K2014 Pin L&H
107248	Motor 10K2014 Pin R&H






SUNRISE MEDICAL
201-337-0201 FAX 201-337-0202
TYPE: 95C2445008AC265255401
Motor DC 20.0V Brake DC 12V
Serial Range 25 to 1 SPIN 145
Date Code 090521 Rev. C
090501E01-112401 CUBES 15.405



Label with Details including Brush Size

LEAD WIRE CHART		
COLOR CODE	DESCRIPTION	FUNC.
BRAKE		
WHITE	BRAKE CONNECTOR	0300
WHITE	BRAKE SWITCH	0300
MOTOR		
RED	+	+
BLACK	-	-
PIN 107248 = 12V OUTPUT BRUSH ROTATION		



CONFIDENTIAL AND PROPRIETARY
 Duplication or Distribution Prohibited
 6/20/12

55

S-636: New Motors

Small Lever Release
New Motors
Geared Freewheel



New Motors feature easy to release levers



New motors are easy to maintain and do not leak grease



CONFIDENTIAL AND PROPRIETARY
 Duplication or Distribution Prohibited
 6/20/12

56

S-636: Dynamic Suspension






CONFIDENTIAL AND PROPRIETARY
 Duplication or Distribution Prohibited
 6/20/12

57

S-636: Full Suspension Base



- The Quickie S-636 offers independent suspension at every point of contact with the ground. This virtually guarantees a smooth ride over rough terrain.



CONFIDENTIAL AND PROPRIETARY
Duplication or Distribution Prohibited
6/25/12
58

S-636: Traction Through Transitions

- The dynamic anti-tip suspension can move upward to keep the rear wheels in contact with the ground during transitions.



CONFIDENTIAL AND PROPRIETARY
Duplication or Distribution Prohibited
6/25/12
59

Key Features Example: Quickie QM-710


- 1 Full Featured Complex Rehab Power Chair**
 - Modular Components Available in 3 Models
 - Industry Proven Components
 - Sleek, Contoured Design
- 2 High Performance Base**
 - SpiderTrac Suspension
 - Traction and Articulation
 - Mid-wheel Drive Intuitive Handling
 - Designed with Minimal Parts Complexity
- 3 Comfortable and Modular Seating**
 - Parts Commonality with Existing Group 3 Chairs
- 4 Versatile Electronics (What's New on the QM-710)**
 - Non-expandable Electronic Option for Challenging Funding
 - Optimized Wiring
 - Assignable Buttons (R-Net Only)




CONFIDENTIAL AND PROPRIETARY
Duplication or Distribution Prohibited
6/25/12
60

QM-710 – Battery Access



- The Batteries connect up in series
- The battery in the front has a single plug Red
- The back battery has the series connection and the Black connector for control wiring harness connection.




Front Battery



Rear Battery



CONFIDENTIAL AND PROPRIETARY
Duplication or Distribution Prohibited
6/20/12

61

QM-710 - Full Featured Complex Rehab Power Chair

Modular Components Available in 3 Models

- Quickie QM-710 (BC) – No Power Option Base with ASAP II Rehab Seat
- Quickie QM-710 (SC) – Single Power Option Tilt (SPOT) Base with ASAP II Rehab Seat
- Quickie QM-710 (MPC) – Multiple Power Option Base (SPOT Tilt and Recline) with Recline Rehab Seat










CONFIDENTIAL AND PROPRIETARY
Duplication or Distribution Prohibited
6/20/12

62

QM-710: Motors

- **Industry Proven Components - Linux Motors**
- Less than 1% return rate, lowest return rate in company history
- Same proven motor as used on the S-636, Pulse 6, and Rhythm
- Use unique dual end-of-line inspection system: both at vendor and at Quickie factory



CONFIDENTIAL AND PROPRIETARY
Duplication or Distribution Prohibited
6/20/12

63

QM-710: Sleek, Contoured Design

- Shroud Available in 5 Colors





CONFIDENTIAL AND PROPRIETARY
Duplication or Distribution Prohibited
8/25/12

64

QM-710: Spider Track Suspension

- Stable, smooth transitions across multiple terrains
 - Each wheel independently suspends






CONFIDENTIAL AND PROPRIETARY
Duplication or Distribution Prohibited
8/25/12

65


QM-710: Spider Track Suspension


- Gas springs provide a dampen to reduce the jolt of an impact
- Example of a spring suspension
- Example of a dampened suspension

Pure Spring Suspension
Very Bouncy



Classic Cadillac Dampened Suspension
Very Smooth





CONFIDENTIAL AND PROPRIETARY
Duplication or Distribution Prohibited
8/25/12

66

22

QM-710: Spider Track Suspension

- When you slowly lower in a pneumatic office chair, it is the gas release that dampens the fall and/or impact.

Pneumatic Office Chair – classic example of a dampening system



CONFIDENTIAL AND PROPRIETARY
Duplication or Distribution Prohibited
6/25/12
67

Triggered Gas Spring Dampening

- PREVENTS collapse of the front casters
- PREVENTS shock associated with non-suspension fixed front casters



CONFIDENTIAL AND PROPRIETARY
Duplication or Distribution Prohibited
6/25/12
68

Triggered Gas Spring Dampening


- Gas dampening system, prevents collapse but dampens the impact
- User maintains seated position during transitions
 - Dampened contact reduces tonal responses
 - Helps prevent user displacement
 - Helps ensure a stable seating posture




CONFIDENTIAL AND PROPRIETARY
Duplication or Distribution Prohibited
6/25/12
69

QM-710: Traction & Articulation

- Climbs 3" curbs with a 22° transition angle
 - Design maximizes front caster articulation range







CONFIDENTIAL AND PROPRIETARY
 Duplication or Distribution Prohibited
 6/20/17

70

QM-710: Traction & Articulation

- Linked suspension moves casters and drive wheels in opposite directions to minimize high centering





CONFIDENTIAL AND PROPRIETARY
 Duplication or Distribution Prohibited
 6/20/17

71

QM-710: Designed With Minimal Parts Complexity

- 33% reduction in parts compared to legacy Quickie high-end power chairs
- Fewer parts reduces
 - Fewer parts to potentially fall out of tolerance, which reduces potential breakdowns
 - No cable adjustments
 - Fewer parts to service creates a reduced service turnaround
 - Easier for technicians to understand






CONFIDENTIAL AND PROPRIETARY
 Duplication or Distribution Prohibited
 6/20/17


72

SpiderTrac Suspension: Gas Strut

- Suspension mechanism accomplishes two tasks with one mechanism
- Suspension prevents forward pitching and high centering with a single mechanism



The linkage between the front caster and rear caster prevent high-centering. The same gas spring is used to absorb shock for both casters.



The rear casters rotate downward as the chair goes over a curb, this triggers a constriction of gas flow in the gas spring, which dampens front caster impact and prevents caster collapse.

STEPS
Science Training & Education Programs

CONFIDENTIAL AND PROPRIETARY
Duplication or Distribution Prohibited
6/25/12




73

Removal of Pneumatic Gas Piston

Gas Spring Removal

- 1. Remove the wheel assembly by using a 5/8" deep socket wrench to remove the 3 bolts that hold the wheel and rim in place (figure 1). Before removing the Gas Spring, raise the center of the chair with a lift so that the suspension can come all of the way down.
- 2. For Gas Spring Removal, Remove the Lower holding bolt (A), using a 5mm Hex key (figure 2)
- 3. Remove the top holding Bolt with a 6mm hex in Pic. 3

Note: The shocks used in this power base are pre-adjusted at the factory.




STEPS
Science Training & Education Programs

CONFIDENTIAL AND PROPRIETARY
Duplication or Distribution Prohibited
6/25/12

74

Removal of Pneumatic Gas Piston

- 4. Once the top holding bolt is removed, loosen the set screw (A) as shown by using a 4mm Hex wrench (figure 4)
- 5. Once setscrew is loosened, the gas spring shaft can be unscrewed and removed from the securing point as shown in (figure 5)
- 6. To re-install the Gas Spring, Reverse the previous procedures. Be sure to adjust the the Gas Springs threaded shaft to the correct distance of 12-13mm for required installation clearance (C) as shown in (figure 6).
- Then tighten the Jam Nut (D) to secure the shaft at the correct distance, and replace the gas spring insert.

STEPS
Science Training & Education Programs

CONFIDENTIAL AND PROPRIETARY
Duplication or Distribution Prohibited
6/25/12


75

Quickie QM-715 HD and QM-720

Quickie QM-715 – 400 lbs
Heavy Duty Package – 5 mph



Quickie QM-720 – 300 lbs
High Speed Package



116515	WTRGBX 4 POLE 6 MPH QM710 LH	Left QM-710
116516	WTRGBX 4 POLE 6 MPH QM710 RH	Right QM-710
120636	WTRGBX 60PH 4POLE QM715 LH	Left QM-715 or QM-710 High Torque
120637	WTRGBX 60PH 4POLE QM715 RH	Right QM-715 or QM-710 High Torque
120638	WTRGBX 130PH 4POLE QM720 LH	Left QM-720 or QM-710 Community Package
120639	WTRGBX 130PH 4POLE QM720 RH	Right QM-720 or QM-710 Community Package




Staircase Training & Education Program


CONFIDENTIAL AND PROPRIETARY
Duplication or Distribution Prohibited
6/20/12

76

Battery Charger – BATC8-CT

- Operation of Charger
 1. Green Flash - indicates power is on
 2. Orange Flash - indicates pre-charge
 3. Orange – indicates charging
 4. Green & Orange Flash – 80% charge rate
 5. Green – fully charged






BATTERY CHARGER INSTRUCTION

1. APPEARANCE

1. Power Cord
2. Charger Plug to Battery
3. Indicator
4. Battery Flash - Power On
5. Orange Flash - Pre-Charge
6. Orange - Charging
7. Green & Orange Flash - Charge 80%
8. Green - Full Charge
9. Red Flash - Fault





Staircase Training & Education Program


CONFIDENTIAL AND PROPRIETARY
Duplication or Distribution Prohibited
6/20/12

77

Back of Battery Charger







Staircase Training & Education Program

CONFIDENTIAL AND PROPRIETARY
Duplication or Distribution Prohibited
6/20/12

78

R-net Expandable Electronics

- Provides simple programming options and outstanding reliability
- Features color screens on the joystick and OMNI display with a consumer clock
- Allows completely unique programming of profiles including the ability to mix proportional and non-proportional controls on one chair




VR-2 Non-expandable Electronics

- Penny & Giles VR-2
 - Programmable
 - With hand held or PC programmer
 - 1 profile with 5 speed increments
 - 5 profiles with 1 programmed speed
 - 70 amp controller
 - Uses the same parameters as Qtronix
 - Provides activation of 2 power options through the joystick
 - Controller to be used with a joystick only




Pulse 6: VR2 (PG Drives) System



VR-70 Controller

- Connector View



The diagram shows the rear of the VR-70 Controller with various connectors labeled. A 'Bus Connector' is at the top left. Below it are 'NH2', 'A1', and 'A2' connectors. To the right is an 'OBC' connector. Further right are 'Motor Connectors' labeled 'M1' and 'M2', with '+' and '-' polarity markings. At the bottom right is a 'Power Connector'.

Bus Connector

Motor Connectors

Power Connector

NH2

A1

A2



OBC

M1


M2

+

-


[illegible]

Tie Wrap Wiring for VR-2




Outside rear view
VR-2 shown

This diagram shows the external wiring of the VR-2 unit. Labels with leader lines identify the following components: Cable Tie (top left), Motor Harnesses (top center), Controller Power Harness (top right), Cable Tie (top right), and Blue Harness (bottom left).



Inside rear view
VR-2 shown



This diagram shows the internal wiring of the VR-2 unit. Labels with leader lines identify the following components: Blue Harness (top left), Controller Power Harness (center left), Motor Harnesses (center right), Cable Tie (top right), Cable Tie (bottom left), and Loop Excess Motor Harness Here (bottom center).



STEPS
Source Training & Innovation Programs

CONFIDENTIAL AND PROPRIETARY
Duplication or Distribution Prohibited
#32010

84



R-net EL-90

- Penny & Giles R-net EL
 - Fully programmable, 5 profiles
 - With DTT or PC programmer
 - 90 amp controller
 - Uses the same parameters as R-net
 - Standard LED joystick (no screen)
 - Compatible with specialty controls




Optional Color Screen Joystick features OBP – On board programming through sequence or dongle for PC usage

Attendant Control : Easily hooks up and has Mode switch usage for other Modes of operation as well




CONFIDENTIAL AND PROPRIETARY
Duplication or Distribution Prohibited
6/20/12
85

R-net EL-90 Control Inputs



Inhibit Line
Off Board Charger Com Connection O.B.C. Motor 1 Power Connect Motor 2



CONFIDENTIAL AND PROPRIETARY
Duplication or Distribution Prohibited
6/20/12
86


R-net Features


- Simple programming with outstanding options
- Color screens
- Unique programming profiles
- Mix proportional and non-proportional controls in one chair.

	FWD	P1	P2	P3	P4
Acc	30	35	40	45	
Acc	20	20			
Dec	40	40			
Dec	30	30			

Global Controls	
Sounder Volume	5
Ending Sleep	n
Act. Entry Axis	1
Chge P1 in Div	yes

OBP MENU	
Controls	>
Latched	>
Actuators	>





CONFIDENTIAL AND PROPRIETARY
Duplication or Distribution Prohibited
6/20/12
87

VR2 Programming

[illegible]

- 1 P = Pilot, P+ = Pilot +
- 2 Not Pilot manufactured before June 1998
- 3 Not Pilot + manufactured before February 2002
- 4 Not V6 manufactured before October 2002
- 5 Not Pilot + manufactured before March 2004
- 6 The V62 parameters are only available on PFI-A versions 9 onwards

VR2 Programming

- Programmer Face



Help Button: Pressing this button displays information regarding the function you have selected. In menus, HELP tells you what each option does. In options, it tells you what to do next.



Up / Yes Button: This button steps up through the menu lists, increases the value of settings and selects functions



Down / No Button: This button steps down through the menu lists, decreases the value of settings and de-selects functions



Enter Button: This button selects options, settings and function states.

VR2 Programming


- **Acceleration**
 - Adjusts the value for forward acceleration of the wheelchair, in increments of 1.
 - The higher the value the faster the acceleration.
- **Deceleration**
 - Adjusts the value for forward deceleration (or braking) of the wheelchair, in increments of 1.
 - The higher the value the faster the deceleration.
- **Turn Acceleration**
 - Adjusts the value for the forward and reverse turn acceleration of the wheelchair, in increments of 1.
 - The higher the value the faster the acceleration.
- **Turn Deceleration**
 - Adjusts the value for the forward and reverse turn deceleration of the wheelchair, in increments of 1.
 - The higher the value the faster the deceleration.

VR2 Programming

- Forward Speed**
 - This sets the MAXIMUM and MINIMUM forward speed of the wheelchair in increments of 1%.
 - There are two available settings.
- Max**
 - The maximum value occurs at the control system's maximum speed setting.
- Min**
 - The minimum value occurs at the control systems minimum speed setting

The value is displayed as a percentage of the wheelchairs total available output.

Therefore if the Max value is set to 80% then the wheelchair will be able to drive at up to 80% of the total available speed when the control system's maximum speed setting is reached.



CONFIDENTIAL AND PROPRIETARY


Duplication or Distribution Prohibited

6/25/12

91

VR2 Programming

- Reverse Speed**
 - This sets the MAXIMUM and MINIMUM reverse speed of the wheelchair, in increments of 1%.
 - There are two available settings.
- Max**
 - The maximum value occurs at the control system's maximum speed setting.
- Min**
 - The minimum value occurs at the control systems minimum speed setting
 - The minimum value is interpreted differently in each control systems.
- Pilot Reverse Speed**
 - This is automatically scaled in relationship to the forward speed setting and calculated as below.
 - Min. reverse = max. forward x min forward
- Pilot+, VSI & VR2 Reverse Speed**
 - This corresponds to the actual value selected with the PP1a Programmer.



CONFIDENTIAL AND PROPRIETARY


Duplication or Distribution Prohibited

6/25/12

92

VR2 Programming

- Turning Speed**
 - This sets the MAXIMUM and MINIMUM turning speed of the wheelchair in increments of 1%.
 - There are two available settings.
- Max**
 - The maximum value occurs at the control system's maximum speed setting.
- Min**
 - The minimum value occurs at the control systems minimum speed setting
 - The value is displayed as a percentage of the wheelchairs total available output.
 - Therefore if the Max value is set to 60% then the wheelchair will be able to drive at up to 60% of the total available speed when the control system's maximum speed setting is reached.
- Power**
 - Reduces the power output of the control system.
 - The parameter is adjustable in steps of 1% between 10% and 100%.
 - The main purpose is to limit damage to furniture or doorways if the wheelchair is being used indoors. This is particularly useful on wheelchairs designed for children. If this parameter is used with the control system programmed for multiple drive profiles, then indoor and outdoor profiles can be set.



CONFIDENTIAL AND PROPRIETARY

Duplication or Distribution Prohibited


6/25/12

93

31

VR2 Programming

- Torque**
 - The Torque parameter boosts the current to the motors at low speed settings. If the motor is stalled, for example, the wheelchair is stuck against an obstacle, such as a door threshold; then this will be automatically detected and the current to the motors will be increased, allowing the obstacle to be overcome.
 - Torque can be set between 0% and 100%
 - A value of 0% means the Torque parameter has no effect. Higher values mean that more current will be permitted in the described stall conditions.
- Tremor Damping**
 - This parameter allows the effects of hand tremor to be reduced. If the user has a condition that results in hand tremor, then increasing the value of Tremor Damping will reduce the effect of the tremor, making the wheelchair more controllable.
 - Tremor Damping can be set between 0% and 100%
 - A value of 0% means Tremor Damping has no effect. Note, even at this value, there is inherent damping in the control system. Higher values apply a higher level of damping.




CONFIDENTIAL AND PROPRIETARY
Duplication or Distribution Prohibited
6/25/12

94

VR2 Programming

- Preset Unit**
 - Selecting this sets all menu parameters to their default values. The default values are stored in the controller by PGDT during manufacturing.
- The Values given prior to the next page are values that are found in the profile of the Program and varies according to the 1st page which gives the different versions of PP1 or QTRONIX programmer versions




CONFIDENTIAL AND PROPRIETARY
Duplication or Distribution Prohibited
6/25/12

95

VR2 Programming

Global Settings

- Sleep Timer**
 - A length of time can be set, such that if the control system accepts no valid input for that period of time, it will power down safely.
 - The time can be set in steps of 1 minute between 0 to 30 minutes.
 - If the value is set to 0, no power down will occur.
- Steer Correct**
 - This factor compensates for any mismatching of motors to ensure that the wheelchair drives directly forward when the control system's joystick is being pushed directly forward.
 - It is normally set to zero but may be varied from -9 to +9 in increments of 1.
 - If the chair is veering to the left, you should increase the setting.
 - If the chair veers to the right, decrease the setting.




CONFIDENTIAL AND PROPRIETARY
Duplication or Distribution Prohibited
6/25/12

96

VR2 Programming

- Joystick Throw**
 - This allows you to program the control system so that full speed can be reached with a reduced joystick movement (throw). This is particularly useful for wheelchair users with limited hand or arm movement.
 - The adjustment can be made manually or by programming actual values.
- Speed Adjust**
 - Sets whether the control system's speed/profile buttons will be active while the wheelchair is driving. You can set this function to Yes or No.
 - Yes Means the control system's speed/profile buttons will be active while driving
 - No Means the buttons will only be active when the joystick is in the central position.
 - This adjustment is particularly useful for users who may accidentally operate these buttons while deflecting the joystick.



CONFIDENTIAL AND PROPRIETARY


Duplication or Distribution Prohibited

6/25/12

97

VR2 Programming

- Invert Joystick**
 - This parameter inverts the direction of travel when moving the joystick.
 - This parameter can be set to On or Off.
 - On Deflecting the joystick Forward will result in Reverse drive.
 - Off Deflecting the joystick Forward will result in Forward drive.
- Actuator Selection**
 - There are two programmable parameters relating to actuator channel 1.
 - Speed**
 - Sets the speed of travel of the actuator connected to channel 1. The speed can be set between 1 and 5. 1 is the slowest, 5 is the fastest.
 - End Force**
 - This value should not normally require any adjustment and should not be altered without the permission of the wheelchair manufacturer. *Amount of current allowed*
 - Actuator 2**
 - The same as one only in reference to actuator channel 2.



CONFIDENTIAL AND PROPRIETARY


Duplication or Distribution Prohibited

6/25/12

98

VR2 Programming

- Actuator End Stop Bleep**
 - This allows the use of an Audible bleep to tell you when the Actuator is stalled at its end stop.
 - This parameter can be set to On or Off
- Speed Adjustment while Driving**
 - This parameter sets whether the speed/profile buttons are active while the wheelchair is being driven. The parameter can be set to on or off.
 - On - Means the buttons are active while the wheelchair is being driven, so the user can make maximum speed setting adjustments (or select a different drive profile) while actually moving.
 - Off - Means the buttons are not active while the wheelchair is being driven, so the joystick must be released and the wheelchair at rest before maximum speed setting adjustments (or different drive profile selections) can be made.



CONFIDENTIAL AND PROPRIETARY

Duplication or Distribution Prohibited


6/25/12

99

33

VR2 Programming

- **Reverse Driving Alarm**
 - Sets whether the control system gives an audible warning while driving in reverse. The parameter can be set to on or off.
 - **On** - Means there is an audible alarm given.
 - **Off** - Means there is not.
- **Profiles**
 - A drive profile is a collection of programmable parameters comprising of Acceleration, Deceleration, Turn Acceleration, Turn Deceleration, Forward Speed, Reverse Speed, Turn Speed and Power. The number of drive profiles is determined by this parameter.
 - If the value is 0, there is one setting for each of the parameters listed previously, and the control system's maximum speed setting can be changed on the normal way.
 - If the value is 2 to 5, there is a corresponding number of drive profiles and each listed parameter can be individually set within a profile. The normal method of maximum speed adjustment can then be used to switch between the available profiles.




Senior Training & Education Programs

CONFIDENTIAL AND PROPRIETARY
Duplication or Distribution Prohibited
6/20/12

100

VR2 Programming

- **Read System Log**
 - The control systems have a diagnostics log facility which stores the number of occurrences of the last eight trip codes.
 - This allows you to view the contents. The display format is as below.
 - 1: Code 2C00, #1
 - 2: Code 3C00, #3
 - No more entries
 - This reads line by line as.
 - Line 1 - trip code 2C00 has occurred once
 - Line 2 - trip code 3C00 has occurred three times
 - Only two trip types recorded.
- **Read Timer**
 - Control has a timer which records how long the wheelchair is in use.
 - The timer runs whenever the joystick is moved away from the center position, and stops when the joystick is returned. The timer records the number of hours the wheelchair has been in use.



Senior Training & Education Programs

CONFIDENTIAL AND PROPRIETARY
Duplication or Distribution Prohibited
6/20/12

101

R-net Expandable Electronics



Joystick Modules



Attendant Module



Power Modules



Programming Key



Specialty Controls



Bus Expansion



Output Modules



Joystick Options



Handheld Programmer (DTT)



PC Programming



Senior Training & Education Programs

CONFIDENTIAL AND PROPRIETARY
Duplication or Distribution Prohibited
6/20/12

102

Diagnostic Test Tool

- The Diagnostic Test tool has been designed specifically to assist technicians, engineers and specialists. The DTT allows users to access, program and share PG Drive Technology controller specific files in a convenient and effective way.
- The Diagnostic Test Tool will allow users to:
 - Program PG Drives Technology controllers
 - Read and save program files
 - Read and save diagnostic logs
 - Read real-time controller related information
 - Save and delete program files locally
 - Save program files to external storage devices, such as memory sticks
 - Manage program + system log files from PC's. When connected the DTT appears as a drive on the PC and all normal PC functionality exists.




Diagnostic Test Tool

- The DTT (Diagnostics Test Tool) has USB port. Files can easily be transferred using a USB drive.
- Holds up to 1000 profiles



DTT Connections


- Two Cables
 - Neutrick Connector for VR2
 - R-net Connector for RNet
- USB an memory stick allows for file transfers and back-up of chair programming



R-net Programming

- On-Board & PC Programming

- On-Board Programming (OBP) uses the graphics LCD on the R-net Joystick Module to provide a clear display of multiple profiles, and allows for easy navigation and adjustment of the screens, using intuitive combinations of joystick movements and button presses. Access to OBP can be made secure via a keycode or a hardware key (dongle).



PC Programming
& OBP through Dongle



CONFIDENTIAL AND PROPRIETARY
Duplication or Distribution Prohibited
8/20/12
106

R-net On-board Programming

- Mode Key Sequence



- If the keycode method is set, then the following button sequence will allow entry to OBP:
 - Note, a dongle will also allow OBP access if this method is set.
- 1. Hold down the Horn button and then hold down the On/Off button until there is a short bleep.
- 2. Power-up sound will occur prior to this bleep.
- 3. Release the Horn button, but continue to hold down the On/Off button until there is a further short bleep.
- 4. Release the On/Off button, there will now be a longer bleep and OBP mode will be entered.



CONFIDENTIAL AND PROPRIETARY
Duplication or Distribution Prohibited
8/20/12
107

R-net OBP Through Mode Switch


- If the hardware dongle method is set, then a PGDT OBP Dongle needs to be connected to allow access to OBP. To use the dongle, follow the sequence below:
 - Turn off the control system.
 - Insert the R-net Dongle along the communication cables in the system configuration. See figure on Right side of Page
 - Turn on the control system
 - After initialization press the Mode key until the OBP screen is reached. See diagram below.



CONFIDENTIAL AND PROPRIETARY
Duplication or Distribution Prohibited
8/20/12
108

R-net OBP Programmable Speeds

FORWARD	
For Accel Max	Sets Maximum Forward Acceleration
For Accel Min	Sets Minimum Forward Acceleration
For Decel Min	Sets Maximum Forward Deceleration
For Decel Max	Sets Minimum Forward Deceleration
For Speed Max	Sets Maximum Forward Speed
For Speed Min	Sets Minimum Forward Speed
REVERSE	
Rev Accel Max	Sets Maximum Reverse Acceleration
Rev Accel Min	Sets Minimum Reverse Acceleration
Rev Decel Min	Sets Maximum Reverse Deceleration
Rev Decel Max	Sets Minimum Reverse Deceleration
Rev Speed Max	Sets Maximum Reverse Speed
Rev Speed Min	Sets Minimum Reverse Speed




CONFIDENTIAL AND PROPRIETARY
Duplication or Distribution Prohibited
©2010

109

R-net OBP Programmable Speeds

TURN	
Turn Accel Max	Sets Maximum Turning Acceleration
Turn Accel Min	Sets Minimum Turning Acceleration
Turn Decel Min	Sets Maximum Turning Deceleration
Turn Decel Max	Sets Minimum Turning Deceleration
Turn Speed Max	Sets Maximum Turning Speed
Turn Speed Min	Sets Minimum Turning Speed
DRIVE	
Power	Reduces power to minimize risk of interior fittings damage
Torque	Torque boost to overcome obstacles at low speed settings
Tremor Damping	Adjustable damping to reduce the effect of hand tremor




CONFIDENTIAL AND PROPRIETARY
Duplication or Distribution Prohibited
©2010

110

R-net Programming

GLOBAL CONTROLS	
Steer Correct	Adjusts the PM output's to compensate for mis-matched motors
Sound Volume	Sets volume of audible feedback from JSM
Endstop Bleep	Sets whether there is a bleep when a seat sub reaches endstop
Act. Entry Axis	Sets the default axis when seat control mode is entered
Change PFI in Drv	Sets whether profile changes are permissible while driving
Speed Adjust	Sets whether the speed buttons on the JSM are active
Speed Adj in Drv	Sets whether speed setting changes are permissible while driving
Momentary Screens	Sets Whether large screens appear at profile and speed changes
Rev Driving Alarm	Sets if the reverse driving alarm is active
EM Stop Switch	Allows latched operation without the use of an EM Stop Switch
Lock Function Enabled	Sets how the lock function is activated
Display Speed	Sets how the speedometer is displayed (miles or kilometers per hour)
Max Display Speed	Sets the operation of the graphical speed display
Power - Up Mode	Sets the Mode that will be active when the system is powered up



CONFIDENTIAL AND PROPRIETARY
Duplication or Distribution Prohibited
©2010

111

R-net Programming

PROFLED CONTROLS

Sleep Timer	Sets the time of inactivity before the system goes to sleep
Chge Mode in Drv	Sets whether mode changes are permissible while driving
Background	Sets the default background for each profile joystick

JOYSTICK

Active Throw	Sets joystick throw via joystick movements
Throw Detail	Sets joystick throw via programming
Active Orientation	Sets joystick orientation via joystick movements
Orientation Detail	Sets joystick orientation via programming
Active Orientation	Sets joystick orientation via joystick movements
Orientation Detail	Sets joystick orientation via programming
Deadband	Sets the joystick deadband (size of neutral position)



CONFIDENTIAL AND PROPRIETARY
Duplication or Distribution Prohibited
8/20/12
112

R-net Programming

STANDBY

Switch to Standby	Sets whether an external button can be used to enter Standby Mode
Standby Time	Sets the time of inactivity before Standby Mode is entered
Mode Select	Sets whether other modes can be selected from Standby Mode
Remote Select	Sets whether a profile can be selected from Standby Mode

LATCHED

Drive	Selects latched drive operation
Actuators	Selects latched actuator operation
Timeout	Sets the timeout period for latched operation
Timeout Sleep	Sets whether a sleep occurs as the timeout period approaches



CONFIDENTIAL AND PROPRIETARY
Duplication or Distribution Prohibited
8/20/12
113

R-net Programming

SM ACTUATORS

Away Speed	Sets actuator axis speed in the up direction
Home Speed	Sets actuator axis speed in the down direction

ISM ACTUATORS

Away Speed	Sets actuator speed in the up direction
Home Speed	Sets actuator speed in the down direction
Acceleration	Sets actuator acceleration
Deceleration	Sets actuator deceleration



CONFIDENTIAL AND PROPRIETARY
Duplication or Distribution Prohibited
8/20/12
114

DTT Programming

File Management




Programming Menus




CONFIDENTIAL AND PROPRIETARY
Duplication or Distribution Prohibited
6/30/12
115

DTT Diagnostics

Immediate Diagnosis



System Tests



System Logs



CONFIDENTIAL AND PROPRIETARY
Duplication or Distribution Prohibited
6/30/12
116


Don't Forget!

- Please fill out evaluation
- Please take Test for CEU's
- Please fill in SSN's
- Have a great day!

**Quickie Pulse 6 –
QM-710 – and S636**
... simple, reliable



CONFIDENTIAL AND PROPRIETARY
Duplication or Distribution Prohibited
6/30/12
117











Sunrise Training & Education Programs

Thank You For Attending!


Ken Kalinowski, Senior Service Technician
Ken.kalinowski@sunred.com

MANUALPOWERADULTPEDIATRICSEATINGGERIATRICSCONTROLSFUNDING

Appendix A


Specialty Control Parts




115207

DTT Kit Includes Cables

VR2 and R-Net



115211 Harness DTT - Neutik



111423 Cable DTT R-Net

STEPS

CONFIDENTIAL AND PROPRIETARY
Duplication or Distribution Prohibited
8/22/12

119

Appendix A

S636


DESCRIPTION	QTY	UNIT PRICE	TOTAL PRICE
115207 DTT KIT INCLUDES CABLES	1	115.00	115.00
115211 HARNESS DTT - NEUTIK	1	115.00	115.00
111423 CABLE DTT R-NET	1	115.00	115.00
SUBTOTAL			345.00
TAXES & FEES			0.00
TOTAL DUE			345.00

FOR OPTIONAL HAND CONTROLS, SEE: PERMETRON/REMOTEDS AND POWER

115207 DTT KIT INCLUDES CABLES

115211 HARNESS DTT - NEUTIK

111423 CABLE DTT R-NET



115207 DTT KIT INCLUDES CABLES

115211 HARNESS DTT - NEUTIK

111423 CABLE DTT R-NET

STEPS

CONFIDENTIAL AND PROPRIETARY
Duplication or Distribution Prohibited
8/22/12

120

Qtronix Joystick Change


Functionality:
The new Pilot + joystick only is configured in two configurations with "1/8" Jack Ports" (118789) or with "Speed Pot & On/Off Switch" (118788).

The new Pilot + joystick mounts differently than the old Qtronix joystick and because of this the mounting plate needs to be changed on Fixed and Center mount Joystick Mounts and the entire Joystick arm needs to be replaced when a swing away mount is utilized. Please see below for replacement components depending on the joystick arm style.

Does your joystick utilize a Round or a Square tube?

Round Tube.....

Fix Mount you need to order		
481211	Q-tronix to Pilot+ adapter Bracket	Qty - 1
220120	8-32 x 1/2 Socket Head Screw	Qty - 2
Right hand Swing Away Mount you need to order		
882498	Complete Swing Away Joystick Mount	Qty - 1
220120	8-32 x 1/2 Socket Head Screw	Qty - 2
Left hand Swing Away Mount you need to order		
882499	Complete Swing Away Joystick Mount	Qty - 1
220120	8-32 x 1/2 Socket Head Screw	Qty - 2
Center Mount you need to order		
481489	Q-tronix to Pilot+ adapter Bracket	Qty - 1
220120	8-32 x 1/2 Socket Head Screw	Qty - 2
221090	1/4-28 x 3/8 Button Head Screw	Qty - 2



CONFIDENTIAL AND PROPRIETARY
Duplication or Distribution Prohibited
8/28/12

124

Qtronix Joystick Change

Square Tube.....


Fix Mount you need to order		
481211	Q-tronix to Pilot+ adapter Bracket	Qty - 1
220120	8-32 x 1/2 Socket Head Screw	Qty - 2
Right hand Swing Away Mount you need to order		
118824	Complete Swing Away Joystick Mount	Qty - 1
220120	8-32 x 1/2 Socket Head Screw	Qty - 2
Left hand Swing Away Mount you need to order		
118825	Complete Swing Away Joystick Mount	Qty - 1
220120	8-32 x 1/2 Socket Head Screw	Qty - 2
Center Mount you need to order		
481489	Q-tronix to Pilot+ adapter Bracket	Qty - 1
220120	8-32 x 1/2 Socket Head Screw	Qty - 2
221090	1/4-28 x 3/8 Button Head Screw	Qty - 2

Qtronix PN	Description	Pilot PN	Description	Feature Comparison				
				Pot	Toggle	Jackie	Right Hand	Left Hand
139515	REMOTE BOX W/1/8 JACK & SS25	118788	JOYSTICK IN POT/TOG NO JACKS	-	-	NA	-	NA
139516	REMOTE BOX W/1/8 JACK & SS25	118789	JOYSTICK IN JACKING POT/TOG	-	-	-	-	NA
139517	REMOTE BOX W/1/8 JACK & SS25	118788	JOYSTICK IN POT/TOG NO JACKS	-	-	NA	NA	-
139518	REMOTE BOX W/1/8 JACK & SS25	118789	JOYSTICK IN JACKING POT/TOG	-	-	-	NA	-

Feature Retained

Feature lost

Feature not available on Qtronix JS



CONFIDENTIAL AND PROPRIETARY
Duplication or Distribution Prohibited
8/28/12

125
